Test	Qualitative assay for the detection of Middle East Respiratory Syndrome Coronavirus
Description	(MERS-CoV) nucleic acid in clinical specimens
Test Use	To aid in the diagnosis of Middle East Respiratory Syndrome in symptomatic persons
Test	Virology
Department	Phone: (860) 920-6662 FAX: (860) 920-6661
Methodology	Real-time Reverse Transcriptase Polymerase Chain Reaction (rRT-PCR)
Availability	Monday-Friday
Specimen	To increase the likelihood of detecting infection, collection of all three specimen types
Requirements	is recommended
	<ul> <li>Lower respiratory specimens (preferred), such as bronchoalveolar lavage,</li> </ul>
	tracheal aspirate, pleural fluid, or sputum
	<ul> <li>Nasopharyngeal and oropharyngeal swabs or washes</li> </ul>
	Serum
Collection	Category B shipping box with cold pack
Kit/Container	To obtain collection kit, refer to Collection Kit Ordering Information
	<ul> <li>Collect lower respiratory specimens within 7 days of symptom onset and before antiviral medications are used. Submit 2-3 mL in a sterile screw capped</li> </ul>
	container.
Collection	<ul> <li>Nasopharyngeal/oropharyngeal swabs submitted in 2-3 mL viral transport</li> </ul>
Instructions	media. Use only synthetic fiber swabs with plastic shafts. Do not use calcium
	alginate swabs or swabs with wooden shafts. NP/OP specimens can be
	combined, placing both swabs in the same vial.
	<ul> <li>Nasopharyngeal or nasal aspirates / washes: collect 2-3 mL into a sterile screw</li> </ul>
	capped collection container.
	<ul> <li>Collect blood by standard venipuncture in a serum separator tube. Submit ≥0.5</li> </ul>
	mL serum, optimally during the first 10-12 days after symptom onset
Specimen	Refrigerate specimens at 2-8° C and transport with an ice pack coolant within 3 days of
Handling &	collection.
Transport	If exceeding 72 hours, freeze at -70°C and ship on dry ice.
Unacceptable	Unlabeled specimen
Conditions	Specimens that have leaked or containers that have broken in transit
	Improperly collected or transported specimens
Requisition	Clinical test requisition, select MERS CoV (Novel Coronavirus)
Form	
	Name and address of submitter (and/or Horizon profile #)
Required	Patient name or identifier, date of birth, town of residence (city, state, zip)
Information	Specimen type or source of collection, test requested
	Please ensure patient name on requisition matches that on the specimen.
Limitations	Testing is limited to those patients symptomatic of respiratory infection and who have
	traveled from the Arabian peninsula or neighboring countries 14 days before illness
	onset. Physician must contact the <b>CT Department of Public Health Epidemiology</b>
	<b>Program</b> at 860-509-7994 (Monday-Friday from 8:30 am – 4:30 pm) or 860-509-8000
	(after-hours/weekends) and submit a "MERS Patient Under Investigation (PUI) Form"

Additional
Comments

Serologic testing on a single serum sample collected 14 or more days after symptom onset is currently available at CDC upon request and approval. MERS-CoV serologic test is for research/surveillance purposes and not for diagnostic purposes

Revision: 8/30/2018