Synhilic Confirmation	
Syphilis Confirmation	I (VDRL & IP-PA)

	Detection of Treponema pallidum antibodies in human serum or plasma. Test	
Test description	results include a quantitative VDRL titer and a qualitative treponemal-specific	
	assay result.	
	To confirm reactive results of non-treponemal syphilis screening (such as	
Test use	VDRL or RPR); as a diagnostic test in individuals with a nonreactive non-	
	treponemal test result but with symptoms suggestive of late syphilis.	
Test	Virology	
Department	Phone: (860) 920-6662, FAX: (860) 920-6661	
Methodology	VDRL slide flocculation test and Treponema pallidum Particle Agglutination	
	Assay (TP-PA)	
Availability	Daily, Monday-Friday	
Specimen	1-2 mL serum (preferred), or	
requirements	plasma collected using EDTA, Sodium Citrate or Heparin anticoagulants	
Collection	To obtain a collection kit, refer to Collection Kit Ordering Information.	
Kit/Container		
Collection	Standard venipuncture technique	
instructions		
Specimen	Store specimen at 2-8° C. Specimens must be received within 5 days of	
Handling &	collection.	
Transport	Transport with ice pack coolant (preferred) or at ambient temperature. Avoid	
	temperature extremes.	
Unacceptable	Unlabeled specimen	
Conditions	Specimen that has leaked or container that has broken in transit	
	Serum that is hemolyzed or chylous	
Requisition	Clinical Test Requisition (select Syphilis Confirmation (VDRL & TP-PA))	
Form		
	Name and address of submitter (and/or Horizon profile #)	
Required	Patient name or identifier, date of birth, town of residence (city, state, zip)	
Information	Specimen type or source of collection, date of collection, test requested	
	Please ensure patient name on requisition matches that on the specimen.	
	False positive TP-PA reactions occur in a small percentage (<1%) of normal or	
Limitations	healthy individuals. False positives may occur in patients with other	
	underlying conditions.	
Additional	Test results include both VDRL and TP-PA. All VDRL results are titered to	
Comments	endpoint.	
Devision deter 0/25		

Revision date: 8/25/15